The logo for CAT, consisting of the letters 'CAT' in a bold, yellow, sans-serif font on a black rectangular background.

Ingenieurbüro CAT  
M. Zipperer GmbH

## Instruction Manual

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---

**Hotplate**  
**Heizplatte**  
**H 30/30D**

# Table of Contents

1	User Instructions .....	4
1.1	Important Instructions for your safety .....	4
2	Danger symbols in this operating manual .....	4
2.1	Danger symbols:.....	4
2.2	Danger levels.....	5
2.3	Warning signs on the device .....	5
3	General safety warnings and instructions .....	5
4	General Information.....	7
5	Scope of Delivery.....	7
6	Intended Use .....	7
7	Setting up and Starting up the Instrument .....	8
7.1	Unpacking the Instrument .....	8
7.2	Electrical Connection .....	8
7.3	Setting up the Instrument.....	8
7.4	Rear Connections .....	9
7.5	Connection of the Pt100 Probe.....	9
7.6	Probe Stand.....	10
7.7	Extension Connectors .....	10
8	The Control Panel.....	11
8.1	Description of the Control Panel Functions .....	11
9	Heating .....	12
9.1	Temperature Control with external Pt100 Probe .....	13
9.2	Temperature Control without Pt100 Probe (Hot Plate Control) .....	13
9.3	Plate Limit (Heating with Pt100 probe) .....	14
9.4	Safety Temperature .....	14
10	Fahrenheit-Celsius Setting .....	15
11	Timer Controlled Operation .....	15
12	Additional Safety Functions .....	15
12.1	Differential Alarm.....	15
12.2	Out of Liquid Check .....	15
12.3	Failure of External Pt100 Probe .....	16
12.4	Internal Over Heating Monitoring .....	16
13	Error Messages.....	16
13.1	Other Messages .....	17
14	Maintenance and Cleaning .....	17
15	Dismantling, Transport and Storage .....	17
15.1	Dismantling .....	17
15.2	Transport and Storage .....	18
16	Warranty and Liability .....	18
17	Technical Data .....	19
18	Repair Return Form.....	20
34	Konformitätserklärung / Declaration of conformity.....	38

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# 1 User Instructions

## 1.1 Important Instructions for your safety

- Every user must read and understand this manual completely before use. Failure to do so can result in serious injury or death.
- Comply with all safety and accident-prevention regulations applicable to laboratory work.
- Follow general instructions for hazard prevention and general safety instructions, e.g. wear protection clothing, eye protection and gloves.
- This operating manual is part of the product. Thus, it must always be easily accessible.
- This instruction sheet does not purport to address all of the safety problems which might result from the use of this device, chemicals, reagents, apparatus or equipment employed in any specific test or protocols. It is the responsibility of the user to consult their authorized safety advisors and establish appropriate health and safety practices and then determine the application of regulatory limitations prior to use.
- Enclose this operating manual when transferring the device to another place.
- If this manual is lost, please request another one. Please contact your dealer or

Ingenieurbüro CAT  
M. Zipperer GmbH  
Wettelbrunner Str. 6  
D-79282 Ballrechten-Dottingen  
Tel. ++49-7634-5056-800  
Fax ++49-7634-5056-801  
[info@cat-ing.de](mailto:info@cat-ing.de)  
[www.cat-ing.de](http://www.cat-ing.de)

## 2 Danger symbols in this operating manual

The safety instructions in this manual appear with the following danger symbols and danger levels:

### 2.1 Danger symbols:

	<b>Hazard point</b>		<b>Electrical shock</b>
	<b>Risk of fire</b>		<b>Explosion</b>
	<b>Bio hazard</b>		<b>Hot surface</b>
	<b>Chemical hazard</b>		<b>Material damage</b>

## 2.2 Danger levels

	Will lead to severe injuries or death
	May lead to severe injuries or death
	May lead to light to moderate injuries
	May lead to material damage

## 2.3 Warning signs on the device

### **WARNING**

This symbol indicates that it is imperative to read and understand the instruction manual prior to operating the instrument. Please highlight points which require special attention in your field of application so they are not overlooked. Disregard of warnings may result in impairment of serviceability as well as in physical harm to the user.

### **WARNING**

This symbol indicates that there are **hot surfaces** on the device. It refers to the hotplate of the device which can reach up to 500°C during operation causing a hazard resulting in serious burns.

- Avoid contact with the hotplate, vessel and heating medium.
- Avoid storage of flammable substances near the instrument.

## 3 General safety warnings and instructions

	<p> <b>DANGER</b> <b>Risk of explosion.</b></p> <ul style="list-style-type: none"> <li>• Do not operate the device in the vicinity of highly flammable or explosive substances. The instrument is not explosion-proof.</li> <li>• Do not use this device for processing any substances which could generate an explosive atmosphere.</li> <li>• Do not use this device to process any explosive or highly reactive substances.</li> <li>• Do not use flammable or explosive substances near the instrument.</li> </ul>
	<p> <b>DANGER</b> <b>Electric shock as a result of penetration of liquid.</b></p> <ul style="list-style-type: none"> <li>• Do not allow any liquids to penetrate the inside of the housing</li> <li>• Switch off the device and disconnect the power plug before starting cleaning or disinfection work. The On/Off Switch on the device does not disconnect the device from the power source.</li> <li>• Use only closed tubes.</li> <li>• Only plug the device back in if it is completely dry, both inside and outside.</li> </ul>
	<p> <b>WARNING</b> <b>Risk from incorrect supply voltage</b></p> <ul style="list-style-type: none"> <li>• Only connect the device to a AC power source with an protective earth (PE).</li> <li>• Only connect the device to voltage sources which correspondent to the electrical requirements on the type label.</li> </ul>
	<p> <b>WARNING</b> <b>Electric shock due to damage to device or mains cable</b></p> <ul style="list-style-type: none"> <li>• Only connect the device to the mains supply if the device and the mains cable are undamaged</li> <li>• Only use devices that have been properly installed or repaired.</li> <li>• In case of danger, disconnect the device from the mains supply by pulling the power plug from the mains socket or by using the isolating device intended for this purpose (e.g. emergency stop switch)</li> </ul>

	<p><b>⚠️WARNING</b> <b>Lethal voltage inside the device</b></p> <ul style="list-style-type: none"> <li>• Do not open the device.</li> <li>• Ensure that the housing is always closed and undamaged so that no parts inside the housing can be contacted by accident.</li> <li>• The On/Off Switch on the device does not disconnect the device from the power source. Remove the plug from the AC power outlet to disconnect the instrument from the mains supply entirely.</li> <li>• Do not allow any liquids to penetrate the inside of the housing.</li> <li>• Repairs are only to be carried out by trained service technicians.</li> </ul>
	<p><b>⚠️WARNING</b> <b>Damage to health due to corrosive or aggressive chemicals</b></p> <ul style="list-style-type: none"> <li>• Observe all markings on the reagent bottles.</li> <li>• Always check the instrument for leaks and air bubbles. Special attention should be directed to determine that all push-ons, threaded connections and suction tubes are firmly in place before beginning operation. Leaking solutions may endanger persons and materials</li> <li>• Use proper connecting vessels, protective clothing and gloves.</li> <li>• Avoid splashes</li> <li>• When dispensing, maintain a physical distance between the instrument and the body.</li> <li>• Dangerous and fuming chemicals must be dispensed in a fume hood.</li> <li>• Only employ the instrument for the purpose intended by the manufacturer, and particularly within the resistance limits of the instrument. If in doubt, contact your supplier, or the manufacturer's factory representative at the phone number shown at the front page of this operating instruction.</li> <li>• Always use the instrument in such a manner that neither the operator, nor any other person is endangered.</li> </ul>
	<p><b>⚠️WARNING</b> <b>Damages to health due to infectious liquids and pathogenic germs.</b></p> <ul style="list-style-type: none"> <li>• When handling infectious liquids and pathogenic germs, observe the national regulations, the biological security level of your laboratory, the material safety data sheets and the manufacturer's application notes.</li> <li>• Wear personal protective equipment</li> <li>• For comprehensive regulations about handling germs or biological material of the risk group II or higher, please refer to the "Laboratory Biosafety Manual" in its respectively current valid version from the World Health Organisation</li> </ul>
	<p><b>⚠️WARNING</b> <b>Damage to health due to contaminated device and accessories</b></p> <p>In the following cases, sample material can be released:</p> <ul style="list-style-type: none"> <li>- improperly sealed tubes</li> <li>- unstable tubes</li> <li>- high vapour pressure of the content so that the seal of the tubes can spring open</li> <li>- damaged sealings</li> <li>- smashed glass tubes</li> </ul> <ul style="list-style-type: none"> <li>• Only mix in closed tubes</li> <li>• Observe the nationally prescribed safety environment when working with hazardous, toxic and pathogenic samples. Pay particular attention to personal protective equipment (gloves, clothing, goggles, etc.), extraction, and the safety class of the lab.</li> <li>• Decontaminate the device and the accessories before storage and shipping.</li> </ul>
	<p><b>⚠️WARNING</b> <b>Risk of fire</b></p> <ul style="list-style-type: none"> <li>• Do not use this device to process any highly flammable liquids</li> </ul>
	<p><b>⚠️CAUTION</b> <b>Poor safety due to inadequate fixing of the unit</b></p> <ul style="list-style-type: none"> <li>• Ensure that the unit is firmly attached to a solid stand.</li> </ul>

**▲ CAUTION** Poor safety due to incorrect accessories and spare parts.

The use of accessories and spare parts other than recommended by Ingenieurbüro CAT, M. Zipperer GmbH may impair the safety, function and precision of the device. Ingenieurbüro CAT, M. Zipperer GmbH cannot be held liable or accept any liability for damage resulting from the use of incorrect or non-recommended accessories and spare parts, or from the improper use of such equipment.

- Only use accessories and spare parts recommended by Ingenieurbüro CAT, M. Zipperer GmbH

## 4 General Information

Our Hotplates are designed according to protection class 1. They were manufactured and tested according to DIN EN 61010. According to these regulations, the devices are designed to meet the requirements for safe and correct operations. To ensure the proper safety and operational functions of the instrument, the user should follow the instructions and safety guidelines in this manual.

**▲ WARNING**

It is the user's task to find out whether the device is suitable for the application. If in doubt clarify this with your dealer or directly with the manufacturer.

**▲ WARNING**

Please comply with all safety and accident-prevention regulations applicable to laboratory work.

## 5 Scope of Delivery

Unpack the delivered contents carefully and check to see that nothing appears physically damaged or is missing. Your delivery should contain the following:

1 Hotplate H 30/30D (230V Order No.: 60226-0050, 120V Order No.: 60226-0051)

1 Instruction Manual

**Optional:**

1 Pt100 temperature probe (Order No.: 60278-0000)

Stand adapter for H 30/30D (Order No.: 60669-0000)

Stand for H 30/30D (Order No.: 60665-0000)

Probe Clamp (Order No.: 8B00562400)

Clamp (Order No.: 60668-0000)

RS485 cable for connection to PC (Order No.: 61703-0000)

RS232 to USB Adapter (Order No.: 30244-0001)

9 pol. Sub-D extension cord (Order No.: 30275-0051)

## 6 Intended Use

These hotplates are mainly to be used to heat hydrous solutions and oil not containing inflammable parts in glass vessels. They are designed for use in chemical and biological laboratories of industrial enterprises, universities and pharmacies.

To ensure maximum service life, observe the specified ambient conditions (temperature and humidity) and ensure that the instrument is not exposed to a corrosive atmosphere.

The onus is on the user to find out whether the device is suitable for his application. If in doubt clarify this with your dealer or with the manufacturer directly.

## 7 Setting up and Starting up the Instrument

### 7.1 Unpacking the Instrument

Unpack the instrument carefully and check to see that it is not damaged. It is important that any damage incurred during transport be recognized at the time of unpacking. Notify your carrier or forwarding agent immediately in case of such damage.

If the instrument is not damaged and all parts are complete you may start to operate the device after reading the instruction manual.

### 7.2 Electrical Connection

- A EURO AC plug (DIN 49441 CEE 7/VII10/ 16 A 250 V) is standard on all instruments. For North America instruments feature standard US plug 20A 125 V 2 pole 3 wire grounding (NEMA L5-20P).
- When operating the instruments in countries with different AC plug systems use an approved adapter or have a qualified electrician replace the AC plug with an approved model suitable for the country of operation.
- The instrument is earthed as supplied. When replacing the original AC plug, ensure that the earth conductor is connected to the new plug!
- When connecting the instrument to an AC power outlet, ensure that your local supply voltage matches that indicated on the instrument's rating plate.

### 7.3 Setting up the Instrument

Please place the device on a fire-proof horizontal and even surface. A minimum distance of 50cm to inflammable materials should always be observed.

**▲ DANGER** The device may not be operated in explosion-prone areas.

**▲ WARNING** The power cable may not come in contact with the hotplate.


**▲ WARNING** The device is not to be used without supervision.

**▲ WARNING** When connecting the instrument to an AC power outlet ensure that your local supply voltage matches the indication on the instrument.


## 7.4 Rear Connections

## 7.5 Connection of the Pt100 Probe

A DIN Connector on the rear of the instrument is intended for the connection of a Pt100 temperature Probe (Order No. 60278-0000). Pt1000 Probes are not suitable.

Connection of a PT 100 sensor

### **CAUTION**

Ensure that the cable of the Pt100 probe does not come in contact with the hotplate.

### **NOTICE**

Use Pt100 probes with electrically isolated tips only. Probes that allow for an electrical connection to the metal chassis of the device may distort temperature measurement.


### **NOTICE**

Ensure that the Pt100 probe tip is submersed to a depth of at least 50mm in the heating medium to allow accurate measurement.

## 7.6 Probe Stand

On the rear panel you can mount an optional stand adapter (Order No.: 60669-0000), in which you can attach an optional probe stand (Order No. 60665-0000).

*(mounted stand adapter)*

## 7.7 Extension Connectors


These are intended for the connection of external modules or RS485-RS232 adapters for RS485 communication to a Personal Computer. Please contact your supplier or manufacturer for more information and order numbers.

Please note that these connectors are intended for manufacturers accessories and cables only. Connection of other cables (e.g. Network cable) and devices will cause permanent damage to the device and to other connected electronics.

## 8 The Control Panel

The control wheel on the left is used for changing set values. Any values and messages to be displayed are shown on the 4 digit 7-segment display.

(Control Panel of the Hotplate)

### 8.1 Description of the Control Panel Functions

#### **7Segment LED Display**

This LED display is used for showing different values and messages (temperature, timer and error messages)

#### **Standby/On Button**

Switch the device on or off (standby). Switching the unit off all set temperature values will be stored internally.

#### **Button: „Plate On/Off“**

To switch on the heating plate after entering the desired set-values this button needs to be pressed. Switching the plate off is also done by pressing this button.

#### **Button: „Plate Temp“**

Display the actual plate temperature or press again for setting the set point temperature. When used with Pt100 probe: display the actual plate temperature or press again for setting the plate limit temperature.

#### **Button: „Safety Temp“**


Display or set the safety temperature

**Button: „Probe Temp“**  
 Display or set the probe temperature (Pt-100)

**Button: „Timer“**  
 Display or set the remaining timer minutes

**Encoder Wheel with Push Button**  
 This encoder wheel is used to change different set values. Turning clockwise will increment values, while turning anti-clockwise will decrement. Pressing the encoder wheel push button will confirm a set value. Any error messages displayed at start can also be reset by pressing the encoder wheel push button.

**Status LED „Plate On/Off“** This LED is lit when the hotplate is on. It indicates that the unit is trying to reach or keep the entered set temperature for hotplate or probe.

**Status LED „Heater“** This LED indicates that the hotplate is currently being exposed to an energy pulse.

**Status LED „Hot“** When the hotplate is hotter than 60°C this LED lights up.

**Status LED „Timer“** When the timer is activated this LED lights up.

## 9 Heating

The hotplate allows for 2 different modes of operation:

- Heat control with external Pt100 probe temperature control (Recommended)
- Heat control with hotplate only (no Pt100 probe connected)

For exact temperature control it is recommended to use an external Pt100 probe.

	<b>▲WARNING</b> Always switch off the device after use!
	<b>▲CAUTION</b> As long as the Heater LED is flashing or lights, the hotplate is supplied with energy. The Hot LED lights when the hotplate is hotter than 60°C. You may get burned!
	<b>▲CAUTION</b> After switching off the device the hotplate stays hot! Therefore, after switching off the device, a warning message appears on the display if the hotplate exceeds 60°C and the power cord is still plugged in.

## 9.1 Temperature Control with external Pt100 Probe

We recommend the use of an external Pt100 probe (Order No. 60278-0000) to allow for accurate temperature control of the heating medium.

	<b>▲WARNING</b> Please ensure that the cable of the Pt100 probe does not come in contact with the hotplate.
---	---

- **Connect** the Pt100 probe with the DIN connector on the rear of the hotplate unit.
- **Place** the Pt100 probe in the vessel (min. 5cm submersion)

**Press Power On/Standby** to switch on the device. Any error messages displayed at start can be reset by pressing the encoder wheel button.

When „**A**“ (Amount) is shown as the first digit in the display **turn the encoder wheel** until the desired volume (in liters) of the heating medium is shown. Press the encoder wheel button to confirm the volume in the vessel.

**Press Probe Temp** (press again if the **Π** does not blink). Set the desired setpoint temperature value by turning the encoder wheel. Press the *Probe Temp* key again or press the encoder wheel button to confirm the set value or wait 3 seconds.

**Press Plate On/Off** to switch on the hotplate (Status LED above the button will now be lit).

	<b>▲CAUTION</b> If a Pt100 sensor is connected, the set value “Plate Temp” doesn’t refer to the set value of the hotplate, but to the temperature limit of the hotplate (see chapter 9.3-Fehler! Verweisquelle konnte nicht gefunden werden.).
	After disconnecting the Pt100 probe the setpoint of the hotplate temperature will be set automatically to 0°C.

## 9.2 Temperature Control without Pt100 Probe (Hot Plate Control)

For certain applications it could be necessary to control the hotplate temperature instead the probe temperature. This is only possible if no Pt100 probe is connected to the device.

**Press Power On/Standby** to switch on the device. Any error messages displayed at start can be reset by pressing the encoder wheel button.

**Press Plate Temp** (press again if the **H** does not blink). Set the desired setpoint temperature value of the hotplate by turning the encoder wheel. Press the *Plate Temp* key again or press the encoder wheel button or wait 3 seconds to confirm the set value.

**Press Plate On/Off** to switch on the hotplate (Status LED above the button will now be lit).

### 9.3 Plate Limit (Heating with Pt100 probe)

For normal operation, the plate limit temperature should be set to the maximum value to achieve optimal heat-up times. Setting the limit temperature to the maximum value allows the microprocessor to freely select the optimum hotplate temperature to reach the set probe-temperature. However, when operating with flammable liquids it may be required to limit the maximum allowed hotplate temperature to a safe value.

Note: Is the hotplate temperature limited to a too low value the set temperature for the fluid cannot be reached or an extremely long heat-up time will be the case.

If the setpoint of the external Pt100 probe is changed, the plate limit temperature will automatically be set to its maximum value. The new value appears in the display and could be changed to the desired value directly by using the encoder wheel.

To set the hotplate limit temperature when Pt100 probe is connected:

Press **Plate Temp** (press again if the **H** does not blink). Set the desired plate limit temperature by turning the encoder wheel. Press the *Plate Temp* key again or press the encoder wheel button or wait 3 seconds to confirm the set value.

### 9.4 Safety Temperature

A safety temperature can be set to shut down the device if there is an error or to prevent an unexpected heat up of the medium.

The safety temperature is set automatically by the unit when the Pt100 probe or hotplate setpoint temperature is changed. If the unit is operated with an external Pt100 probe and the setpoint is changed the device automatically set the safety temperature 15°C above the new setpoint. If the liquid temperature exceeds the safety temperature (temperature measured with the external probe) the unit will shut down.

Conversely, if the unit is operated without an external probe and the hotplate setpoint is changed the device automatically set the safety temperature 15°C above the new hotplate setpoint. If the temperature measured at the hotplate (internal sensor) exceeds the safety temperature, the unit will also shut down.

If the device automatically changes the safety temperature the new value is shown in the display for a short time. After the device automatically set the new value another safety temperature can be set manually by user.

To set the safety temperature:

Press **Safety Temp** (the **.°** starts to blink). Set the desired Safety temperature by turning the encoder wheel. Press Safety Temp again to confirm the set value or press the encoder wheel button or wait 3 seconds. The set value will be effective until setpoint of *Probe Temp* or setpoint of *Plate Temp* are changed and the *Safety Temp* is automatically set.

When a Pt100 probe is used the safety temperature refers to the temperature measured by the Pt100 probe. Otherwise the safety temperature refers to the hotplate.

In safety critical applications the proper use of the safety temperature is advised. Whenever the setpoint temperature of the external probe (Pt100) or the hotplate is changed the current safety temperature is set to 15°C above the setpoint of the *Probe Temp* or the *Plate Temp* value. To set the safety temp as close to the actual setpoint value as possible, set the probe or plate temperature first, before changing the safety temperature value.

## 10 Fahrenheit-Celsius Setting

This hotplate can show temperature values in degrees Celsius as well as Fahrenheit. When the device is switched on in Celsius mode the normal start sequence will be displayed. In Fahrenheit Mode the message **ΦAHP** will be displayed after switching on the device.


You can change between Celsius and Fahrenheit mode in the following way:

Press **Power On/Standby** to switch on the device.

Shortly **after** switching the device on **press the Probe Temp Key and hold**.

The display now shows a flashing **Y** and you can change between the settings **XEL** (Celsius) or **ΦAH** (Fahrenheit) by turning the encoder wheel. Press the encoder wheel to confirm the selected temperature mode. The unit will remember the selected temperature mode until it is changed again by the user.

## 11 Timer Controlled Operation

This device has an internal timer function. The timer allows for automatic switching off. The switch off time can be selected in the following way.

E.g. a set time of 10 minutes will cause the shut-off of the unit after 10 minutes.

Press **Timer key** (press again if the **X** does not blink). Turn the **encoder wheel** to select the desired switch off time (remaining on-time in minutes). Press the Timer key again or press the encoder wheel button or wait 3 seconds to confirm.

Minimum time until Switch off: 1 minute  
Maximum time until Switch off: 999 minutes (16h 39min)

## 12 Additional Safety Functions

The following safety functions provide additional protection from uncontrolled heating. If a safety problem is detected, the device turns off automatically. The reason for the shutdown is displayed on the LED-display. When switching on the device next time, you have to confirm the displayed safety problem by pressing the encoder wheel button.

### 12.1 Differential Alarm

Breakage of the vessel could create a hazardous situation. To monitor this situation problem this instrument has a *'differential alarm'* feature. The *'differential alarm'* switches the unit off if the probe has fallen out of the liquid (e.g. breakage of glass). The *'differential alarm'* operates if the temperature fall detected by the Pt100 probe exceeds the internal differential alarm sensitivity. The unit then switches off.

### 12.2 Out of Liquid Check

The unit monitors, whether the Pt100 probe tip is immersed in the liquid in the vessel on the hotplate. The microprocessor checks whether the probe temperature changes in relation to the temperature changes of the hotplate. If the microprocessor detects no increase in the liquid temperature over a certain time, although the hotplate temperature is rising, a warning message is shown on the display and the unit shuts down.

#### Note:

This safety feature is only active if the setpoint temperature of the probe is at least 10°C higher than the actual measured temperature.

### 12.3 Failure of External Pt100 Probe

Disconnection or failure of the external Pt100 probe after the device is switched on will be detected and cause the unit to switch off.

### 12.4 Internal Over Heating Monitoring

The temperature inside the device is constantly monitored. A safety mechanism will switch the device off in case the internal temperature exceeds a safe value.

## 13 Error Messages

This unit has an on-board self-diagnostic program that detects possible failures. In case of an error message the program switches off the unit. To enable the user to prevent the problem from re-occurring, the unit will display any of the following error messages as it is switched on again. To reset any error message press the encoder wheel button:

LED Display:	Error name:	Description/Reasons :
HEρP	hotplate failure	The hotplate temperature exceeded its maximum value The hotplate temperature sensor may be broken
ΟΦΕΡ	Off Error	The device was not switched off with the power button. Disconnecting the power cable or a mains power failure may be the source for error.
Η.:ΑΦ	Plate > Safety	The measured hot plate temperature was higher than the set safety temperature for the hotplate (Safety Temp.)
Π.:ΑΦ	Probe > Safety	The measured probe plate temperature was higher than the set safety temperature for the Pt100 probe (Safety Temp.)
ΧΛΟΧ	TIMER expired	Unit switched off after timer expired (→ no error, user information only).
ΔΙΦΦ	Differential Alarm	The differential alarm was triggered The temperature drop in the probe was too quick.
ΠΦΑΙ	Probe failure	The Pt100 probe is broken or damaged The probe was disconnected during operation.
ΤΡΙΑ	Triac failure	The heating plate output stage is damaged.
ΗΦΑΙ	Plate Sensor Fail	The hotplate sensor is broken
ΗΟΥ.:	Housing overtemp	The internal temperature exceeded its maximum safe value
νΛΙ9	Out Of Liquid	The external temperature probe is not inserted into the medium.
ΧΟντ	Contact Thermometer Error	The contact thermometer is damaged or was disconnected during operation.
ΥδτΟ	Watch Dog Timer	The Watch Dog Timer was triggered.

### 13.1 Other Messages

LED Display:	Error name:	Description/Reasons :
Φαηρ	Fahrenheit	°Fahrenheit Mode
Οφφ	OFF	Unit was switched off by ON/OFF key
4850	RS232 Off	Unit was switched off by RS232/RS485 PC Communication
ΗΟτ	Hotplate is hot	The device is switched off but the hotplate exceeds 60°C.

## 14 Maintenance and Cleaning

The outer casing is mostly 304 grade stainless steel, the hotplate consist of aluminium. Together with a chemically resistant splash-proof membrane key-pad and therefore easily cleaned with warm water and any proprietary liquid laboratory detergent. Do not use steel wool or any similar plastic wool sponge or any aggressive cleaning agents to clean the unit.

Crust, limestone and water stains may be cleaned with the usual household detergents. Please remove remaining detergent thoroughly with a damp cloth or sponge, since some of the detergents react caustic with higher temperatures. If soils are already encrusted use a razorblade scratcher. Plastics, kitchen foil as well as stuff containing sugar must be removed immediately.

Mechanical supports are: Razorblade scratchers, non-scratching sponges and non-scratching steel wool.

The H 30/30D should only be opened and repaired by authorized service personnel. Any work on the electronics in the unit should only be carried out by knowledgeable, trained personnel. Any attempt by the user to repair the unit will immediately render the guarantee null and void. Please contact your local distributor in the event of a problem.

	<b>▲WARNING</b> In case of malfunction do not attempt to repair the unit. There are no user-serviceable parts in this instrument.
	<b>▲WARNING</b> Unplug the AC power cord before opening the housing.

## 15 Dismantling, Transport and Storage

### 15.1 Dismantling

1. Switch off the unit by switching the rocker switch (1) to the „OFF“-Position.
2. Set speed to “0” position by turning the speed control knob to “0”.
3. Disconnect the unit from the mains.
4. Now you may remove the instrument from the working area.

## 15.2 Transport and Storage

### Prior to transport:

Switch the instrument off and unplug the power supply.

Place the instrument and its parts in its original packaging or another suitable container to protect it during transport. Close the packaging with adhesive tape.

Store the instrument in a dry environment. Please observe the specified conditions of the ambient (temperature and humidity).

Do not subject the instrument to mechanical shocks or vibration during transporting it.

In case you do not use the original packaging please mark the box with the following notes:

- Glass symbol (handle with care, fragile)
- Umbrella (keep dry)
- Content (list of content)

Please dispose of used instruments and defective components at your local recycling collection point. Prior to disposal, please sort according to materials: metal, glass, plastic, etc. Also be sure to dispose of the packing material in an environmental-friendly manner.

WEEE No. DE 54398924

## 16 Warranty and Liability

The manufacturer agrees to correct for the original user of this product, either by repair, or at the manufacturer's election, by replacement, any defects in material or workmanship which develop within 24 months after delivery of this product to the original user. In the event of replacement, the replacement unit will be warranted for the remainder of the original twenty-four (24) month period of ninety (90) days, whichever is longer.

If this product should require service, contact your distributor or manufacturer for necessary instructions.

This warranty shall not apply if the defect or malfunction was caused by accident, neglect, unreasonable use, improper service, or other causes not arising out of defects in material or workmanship.

There are no warranties, expressed or implied, including, but not limited to, those of merchantability or fitness for a particular purpose, which extended beyond the description and period set forth herein. The manufacturer's sole obligation under this warranty is limited to the repair or replacement of a defective product and the manufacturer shall not, in any event, be liable for any incidental or consequential damage of any kind resulting from use or possession of the product.

### Note

The onus is on the user to find out whether the device is suitable for his application. If in doubt clarify this with your dealer or with the manufacturer directly.

## 17 Technical Data

Type H 30/30D	Specifications
Mains	120/230 Volt; 50/60 Hz (See rear panel)
Display	4 digit 7Segment-LED display
Temperature control	<ul style="list-style-type: none"> <li>- microprocessor controlled (PID Control)</li> <li>- Pt100 probe adjustable in 1°C Steps</li> <li>- Hotplate adjustable in 1°C steps up to 380°C</li> <li>- programmable probe temperature (Pt100)</li> <li>- programmable hotplate temperature</li> <li>- programmable safety temperature</li> </ul>
Safety system	detects and protects against the following hazardous situations: -hotplate failure (hotplate exceeds max. allowed temperature) -probe failure (disconnection or break of Pt100 probe) -detection of out-of-liquid condition, e.g. stirring vessel breaks -differential safety cut-out (external probe senses an extreme drop in temperature of liquid being stirred/heated) -triac failure
Switch-off time delay (Timer)	programmable switch-off time delay min. delay: 1 minute max. delay : 999 minutes
Temperature sensors	Hotplate: Thermocouple External: Platinum Resistance Sensor (Pt100)
Setpoint resolution hot plate	1°C / 1K (1°F in Fahrenheit Mode)
Setpoint resolution Pt100 Probe	1°C / 1K (1°F in Fahrenheit Mode)
Typical temperature stability with Pt100 use*	± 1.0 °C*
Temperature range heating medium	Up to 250°C
Permissible Ambient Temperature	5 - 40 °C
Permissible Humidity	80 % RH
Safety class acc to DIN 40050	IP32
Case dimensions W x D x H	315 x 350 x 145 mm
Weight	7.4 kg
Hot plate material	Aluminium
Hot plate dimensions	300 mm square
Heating power	2000 Watts

\* 2L H<sub>2</sub>O at 60°C Setpoint, Ambient: 23°C. Maximum variance during 60 min test.

## 18 Repair Return Form

<p>Important:</p> <ul style="list-style-type: none"> <li>- Calibration Service is automatically performed for instruments that require adjustments</li> <li>- We contact you by sending a repair estimate.</li> <li>- Repairs will only be carried out <u>after</u> your release.</li> </ul>
--

CONTACT/USER INFORMATION	
Contact:	Phone No.:
Fax No.:	Email:
Billing:	Shipping:
Company	Company
Address	Address

INSTRUMENT INFORMATION	
Model	Serial No.
Please describe all problems/malfunctions	

Operating Conditions (please fill in if applicable)		
Ambient Temp.	Humidity	Speed
Load	Volume	Viscosity
Temperature in °C	Sample Temperature	Operating Time
Sample Description*		

\*NOTE: If the instrument was exposed to hazardous material, it must be decontaminated BEFORE returning it to Ingenieurbüro CAT, M. Zipperer GmbH and an MSDS for hazardous material must be included with the instrument.

RETURN SHIPPING
-----------------

UPS <input type="checkbox"/>	Air Parcel Post <input type="checkbox"/>	Collect** <input type="checkbox"/>	Other** <input type="checkbox"/>
------------------------------	--	------------------------------------	----------------------------------

\*\*Your account number is required for UPS collect respectively the address and contact of your preferred forwarder if you choose any other transport means.

PACKAGING INSTRUCTIONS TO RETURN AN INSTRUMENT FOR REPAIR
---

- Remove all accessories (e.g. homogeniser tools, stirring paddles, stands, clamps, cables) from the instrument
- Clean excess testing material off the instrument/accessory
- Include MSDS sheets for all hazardous materials used with this instrument
- Pack the instrument in its original box. If the box is not available, take care to wrap the instrument and accessories with enough material to support them.
- DO NOT send accessories unless there is a problem with them. When sending back, wrap each item to avoid contact with the instrument.
- Pack the instrument and related items in a strong box for shipping. Mark the outside of the box with handling instructions.

Example: "Handle with care" or "Fragile- Delicate Instrument" and send to

**Ingenieurbüro CAT M. Zipperer GmbH**  
**Service Department**  
**Wettelbrunner Str. 6**  
**D - 79282 Ballrechten-Dottingen**

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We declare under our sole responsibility that this product (see table) corresponds to the regulations and conforms with the standards or standardized documents:

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The agreement with further valid guidelines/regulations following for the product is explained:

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Applied one harmonized standards in particular:

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In the case of a modification of the unit which has not been agreed on with us, this declaration becomes null and void and the warranty expires.

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