The logo for CAT, consisting of the letters 'CAT' in a bold, yellow, sans-serif font on a black rectangular background.


Ingenieurbüro CAT
M. Zipperer GmbH

Operating Instructions

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Roller Mixer / Rollenmischer

RM 5F

Content

1	Scope of Delivery	3
2	User Instructions.....	3
2.1	Important Instructions for your safety	3
2.2	Danger symbols in this operating manual	4
2.2.1	Danger symbols:	4
2.2.2	Danger levels	4
3	General safety warnings and instructions	4
3.1	Warning signs on the device.....	5
4	General Information.....	6
5	Intended Use	6
6	Unpacking.....	6
7	Set-up	7
7.1	Electrical Connection	7
7.2	Operation (Mixing)	7
8	Cleaning and Maintenance.....	8
9	Dismantling, Transport and Storage.....	8
9.1	Dismantling	8
9.2	Transport and Storage.....	8
10	Warranty and Liability	8
11	Technical Data	9
12	Repairs.....	9
12.1	Repair Return Form	10
11	Konformitätserklärung / Declaration of conformity.....	19
12	Bedienungsanleitung / User Instructions Power supply	20

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1 Scope of Delivery

Inspect the content of the package for damage or missing parts:

1 Roller Mixer RM 5F	PN: 60207-0020
1 Instruction Manual	
1 Power Supply	
consisting of:	
Switching Power Supply 24V DC	PN: 40207-0011
Set of input plugs	PN: 30213-0057

2 User Instructions

2.1 Important Instructions for your safety

- Every user must read and understand this manual completely before use. Failure to do so can result in serious injury or death.
- Follow general instructions for hazard prevention and general safety instructions, e.g. wear protection clothing, eye protection and gloves.
- This operating manual is part of the product. Thus, it must always be easily accessible.
- Enclose this operating manual when transferring the device to another place.
- If this manual is lost, please request another one. Please contact your dealer or

Ingenieurbüro CAT
M. Zipperer GmbH
Wettelbrunner Str. 6
D-79282 Ballrechten-Dottingen
Tel.: ++49-(0)7634-5056-800
Fax: ++49-(0)7634-5056-801
www.cat-ing.de
info@cat-ing.de

2.2 Danger symbols in this operating manual

The safety instructions in this manual appear with the following danger symbols and danger levels:

2.2.1 Danger symbols:

	Hazard point		Bio hazard
	Risk of fire		Explosion
	Crushing		Material damage

2.2.2 Danger levels

DANGER	Will lead to severe injuries or death
WARNING	May lead to severe injuries or death
CAUTION	May lead to light to moderate injuries
NOTICE	May lead to material damage

3 General safety warnings and instructions

	<p>DANGER! Risk of explosion.</p> <ul style="list-style-type: none"> Do not operate the device in the vicinity of highly flammable or explosive substances. The instrument is not explosion-proof. Do not use this device for processing any substances which could generate an explosive atmosphere. Do not use this device to process any explosive or highly reactive substances.
	<p>WARNING! Damages to health due to infectious liquids and pathogenic germs.</p> <ul style="list-style-type: none"> When handling infectious liquids and pathogenic germs, observe the national regulations, the biological security level of your laboratory, the material safety data sheets and the manufacturer's application notes. Wear personal protective equipment For comprehensive regulations about handling germs or biological material of the risk group II or higher, please refer to the "Laboratory Biosafety Manual" in its respectively current valid version from the World Health Organisation

	<p>WARNING! Risk of fire</p> <ul style="list-style-type: none"> • Do not use this device to process any highly flammable liquids
	<p>WARNING! Damage to health due to contaminated device and accessories</p> <p>In the following cases, sample material can be released:</p> <ul style="list-style-type: none"> - improperly sealed tubes - unstable tubes - high vapor pressure of the content so that the seal of the tubes can spring open - damaged sealing - smashed glass tubes <ul style="list-style-type: none"> • Only mix in closed tubes • Observe the nationally prescribed safety environment when working with hazardous, toxic and pathogenic samples. Pay particular attention to personal protective equipment (gloves, clothing, goggles, etc.), extraction, and the safety class of the lab. • Decontaminate the device and the accessories before storage and shipping.
	<p>CAUTION! Poor safety due to incorrect accessories and spare parts.</p> <p>The use of accessories and spare parts other than recommended by Ingenieurbüro CAT, M. Zipperer GmbH may impair the safety, function and precision of the device. Ingenieurbüro CAT, M. Zipperer GmbH cannot be held liable or accept any liability for damage resulting from the use of incorrect or non-recommended accessories and spare parts, or from the improper use of such equipment.</p> <ul style="list-style-type: none"> • Only use accessories and spare parts recommended by Ingenieurbüro CAT, M. Zipperer GmbH
	<p>CAUTION! Crush hazard due to moving parts</p> <ul style="list-style-type: none"> • Do not replace any consumables as long as the device is running. • Do not open the coverage as long as the device is running • Only switch-on the device with properly mounted rolls and completely closed coverage • Do not place hands or objects into the area under the rollers while the machine is running.
	<p>NOTICE! Damage to electronic components due to condensation</p> <p>Condensate can form in the device after it has been moved from a cool environment to a warmer environment.</p> <ul style="list-style-type: none"> • After installing the device, wait at last for 3 hours. Only then connect the device to the mains.
	<p>NOTICE! Damage to electronic components due to incorrect supply voltage</p> <p>Incorrect supply voltage may destroy electrical components.</p> <p>Only use a power supply with a 24 V DC, min. 600 mA output</p>
	<p>NOTICE! Damage from the use of aggressive chemicals.</p> <ul style="list-style-type: none"> • Do not use any aggressive chemicals on the device or its accessories such as strong or weak bases, acids, acetone formaldehyde, halogenated hydrocarbons or phenol. • If the device has been contaminated by aggressive chemicals, immediately clean it by means of a mild cleaning agent.

3.1 Warning signs on the device

WARNING!

This symbol indicates to read the instruction manual carefully prior to operation of the instrument. Please mark points which require special attention in your field of application so they are not overlooked. Disregarding of warnings may result in impairment of serviceability as well as impairment of the user.

4 General Information

The roller mixer RM 5 is designed in accordance with Safety Class III (EN 6114), and built and tested in accordance with DIN EN 61010.

According to these regulations, the units are designed to meet the requirements for safe and correct operations. To maintain the proper safety and operational functions of the instrument, the user should follow the instructions and safety guidelines in this manual.

5 Intended Use

The roller mixers are mainly used to mix hydrous solutions not containing inflammable parts in glass or plastic vessels. They are designed for use in chemical and biological laboratories of industrial enterprises, universities and pharmacies.

To ensure maximum service life, observe the specified ambient conditions (temperature and humidity) and ensure that the instrument is not exposed to a corrosive atmosphere.

6 Unpacking

Unpack the instrument carefully and check to see that it is not damaged. It is important that any damage incurred in transport be recognized at the time of unpacking. Notify your carrier or forwarding agent immediately in case of such damage.

▲WARNING

- Read this instruction manual carefully before operating the instrument. Should there be any additional questions, after reading these instructions, concerning the set-up, operation or warranty, please contact either your distributor, or the manufacturer.
- After reading and understanding the instruction manual you may now start operating the unit.
- Store the instruction manual in a place easily accessible to every user.

7 Set-up

⚠ DANGER

- The units are not to be used in rooms with danger of explosion.
- The unit is not to be used without supervision.


Please put the unit on to a fire-proof respectively non-combustible even surface.

7.1 Electrical Connection

The roller mixer RM 5F is powered by a 24 Volt DC power supply with an angled hollow plug 5.5/2.1 mm

The input socket is located on the bottom of the housing, the On/Off switch (1) on the side of the housing.

7.2 Operation (Mixing)

Five synchronised rotating PVC rollers form a homogeneous reel roller bearing for the use of **closed** test tubes. The tubes are both, simultaneously rotated along their longitudinal axis and rocked back and forth on their bearings.

- Place the samples to be mixed between the rollers
- Start mixing function with rocker switch (1)

8 Cleaning and Maintenance

The housing and the 5 PVC-rollers are acid-resistance covered for long lasting use and easy cleaning with warm water and dish washing liquid. Do not use steel wool or a similar scratching sponge to clean the units. The device is maintenance-free.

NOTICE!

- Should repairs become necessary, only trained and authorized service personnel may open and repair the instrument.
- Unplug the power supply before opening the housing.

9 Dismantling, Transport and Storage

9.1 Dismantling

1. Switch off the unit by using switching the rocker switch (1) to the „OFF“-Position.
2. Disconnect the unit from the power supply.
3. Now you may remove the instrument from the working area.

9.2 Transport and Storage

Prior to transport:

- Decontaminate the device and the accessories before storage and shipping.
- Switch the instrument off and unplug the power supply.
- Place the instrument and its parts in its original packaging or another suitable container to protect it during transport. Close the packaging with adhesive tape.
- Store the instrument in a dry environment. Please observe the specified conditions of the ambient (temperature and humidity).
- Do not subject the instrument to mechanical shocks or vibration during transporting it.

In case you do not use the original packaging please mark the box with the following notes:

- Glass symbol (handle with care, fragile)
- Umbrella (keep dry)
- Content (list of content)

Please dispose of used instruments and defective components at your local recycling collection point. Prior to disposal, please sort according to materials: metal, glass, plastic, etc. Also be sure to dispose of the packing material in an environmental-friendly manner.

10 Warranty and Liability

The manufacturer agrees to correct for the original user of this product, either by repair, or at the manufacturer's election, by replacement, any defects in material or workmanship which develop within 24 months after delivery of this product to the original user. In the event of replacement, the replacement unit will be warranted for the remainder of the original twenty-four (24) month period of ninety (90) days, whichever is longer.

If this product should require service, contact your distributor or manufacturer for necessary instructions.

This warranty shall not apply if the defect or malfunction was caused by accident, neglect, unreasonable use, improper service, or other causes not arising out of defects in material or workmanship.

There are no warranties, expressed or implied, including, but not limited to, those of merchantability or fitness for a particular purpose, which extended beyond the description and period set forth herein.

The manufacturer's sole obligation under this warranty is limited to the repair or replacement of a defective product and the manufacturer shall not, in any event, be liable for any incidental or consequential damages of any kind resulting from use or possession of the product.

11 Technical Data

Model	RM 5F
Part. No.	60207-0020
Power Supply	24 V DC, min. 600 mA
Power	15 W
Usable Roll Length (mm)	334
Roll Dia. (mm)	32
Distance from roll to roll (mm)	36
Max. Payload (kg)	5
Speed (RPM)	40
Tilt angle	± 3°
Overall Dimension (WxDxH) (mm)	180 x 480 x 113
Weight (Kg)	3.7

Condition of Surrounding:

Permissible ambient temperature : +5 °C to + 40 °C
Permissible ambient temperature : max 80% relative humidity

12 Repairs

▲WARNING

When returning instruments for repair, which had contact with hazardous substances:

- Decontaminate the device and the accessories before shipping.
- Provide precise information on the relevant medium
- Take protective measures to ensure the safety of our receiving and maintenance personnel
- Mark the package as appropriate for hazardous materials.

12.1 Repair Return Form

CONTACT/USER INFORMATION	
Contact:	Phone No.
Fax No.	Email:
Billing:	Shipping:
Company	Company
Address	Address

INSTRUMENT INFORMATION	
Model	Serial No.
Please describe all problems/malfunctions	

Operating Conditions (please fill in if applicable)		
Ambient Temp.	Humidity	Speed
Load	Volume	Viscosity
Temperature in °C	Sample Temperature	Operating Time
Sample Description*		

*NOTE: If the instrument was exposed to hazardous material, it must be decontaminated BEFORE returning it to Ingenieurbüro CAT, M. Zipperer GmbH and an MSDS for hazardous material must be included with the instrument.

RETURN SHIPPING

UPS <input type="checkbox"/>	Air Parcel Post <input type="checkbox"/>	Collect** <input type="checkbox"/>	Other** <input type="checkbox"/>
------------------------------	--	------------------------------------	----------------------------------

**Your account number is required for UPS collect respectively the address and contact of your preferred forwarder if you choose any other transport means.

PACKAGING INSTRUCTIONS TO RETURN AN INSTRUMENT FOR REPAIR

- ✓ Remove all accessories from the instrument
- ✓ Clean excess testing material off the instrument/accessory
- ✓ Include MSDS sheets for all hazardous materials used with this instrument
- ✓ Pack the instrument in its original box. If the box is not available, take care to wrap the instrument and accessories with enough material to support them.
- ✓ DO NOT send pedestal stand unless there is a problem with the upright rod, clamp or base. If there is a problem with the stand remove the upright rod from the base and individually wrap each item to avoid contact with the instrument. (Applicable for overhead stirrers and homogenisers.)
- ✓ Pack the instrument and related items in a strong box for shipping. Mark the outside of the box with handling instructions.

Example: "Handle with care" or "Fragile- Delicate Instrument" and send to:

Ingenieurbüro CAT, M. Zipperer GmbH
 Service Department
 Wettelbrunner Str. 6
 D-79282 Ballrechten-Dottingen

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We declare under our sole responsibility that this product (see table) corresponds to the regulations and conforms with the standards or standardized documents:

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The agreement with further valid guidelines/regulations following for the product is explained:

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Applied one harmonized standards in particular:

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In the case of a modification of the unit which has not been agreed on with us, this declaration becomes null and void and the warranty expires.

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Power supply

Foreword

This AC/DC adapter is designed as external power supply for electrical equipment.
Please read and follow the instructions before first use.

Installation

Before first use please check the mains voltage to match the input voltage specified on the power supply. Then check the output voltage, its polarity, the output current and power to match the requirements of the equipment to be connected with.

The connection of the plug must be stable; twisting or bending of the cable should be avoided.

Provide adequate ventilation during operation. The adapter should not be covered or used near the heating sources and may not be exposed to direct sunlight.

If your device is not used for a long time we recommend to disconnect the adaptor from main voltage in order to avoid any damage due to overvoltage for example (spikes, lightning).

Safety

Do not use the adapter if it is open or its housing is damaged; furthermore is prohibited to continue using the adapter in case of main connection damage - life danger due to potentially electric shock.

Repairs may be performed only by an authorized technician; any unauthorized modifications or repairs to the unit will void the warranty.

Make sure the DC-connector on secondary side of the adapter is suitable for your device - any forced connection of a non-matching DC connector may damage the device.

Protect the adapter from water and moisture.

Maintenance

For cleaning you may not use cleaning agents. Wipe with a dry cloth. Disconnect the adapter from the mains before cleaning.

Disposal

Electronic devices are hazardous waste and do not belong to the normal house waste. If the device reached end of its life dispose it according to the applicable legal provisions of the municipal collection points. Disposing in the domestic waste is not permitted.

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MEAN WELL ENTERPRISE Co., Ltd.
Schukat electronic Vertriebs GmbH