Instruction Manual

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Shaker / Schüttelgerät S 26

Table of Contents

1	User Instruktion	4
1.1	Important Instructions for your safety.....	4
1.2	Danger symbols in this operating manual.....	5
1.2.1	Danger symbols:.....	5
1.2.2	Danger levels.....	5
1.2.3	Warning signs on the device	5
2	General safety warnings and instructions	6
3	Intended Use.....	8
4	General Information	8
5	Unpacking the Instrument	9
5.1	Scope of delivery	9
6	Set-up and Operation.....	9
6.1	The Front Panel Control	9
6.2	Description of the front panel functions.....	10
6.3	Set-up	10
6.3.1	Electrical Connection	11
7	Operation Shaking.....	11
7.1	Indexing mechanism	11
8	Timer controlled operation.....	12
9	Error Messages.....	12
10	Serial Interface (RS232) of S 26.....	13
10.1	Connection between PC and S 26.....	13
10.2	Format of an RS232 command	13
10.3	RS232 - commands	14
10.3.1	Format of the Controller Handshake	15
10.3.2	RS232-application example	16
11	Maintenance and Cleaning	17
11.1	Cleaning	17
11.2	Maintenance.....	17
12	Dismantling, Transport and Storage	17
12.1	Dismantling.....	17
12.2	Transport and Storage	17

13	Disposal.....	18
14	Warranty and Liability.....	18
15	Repairs.....	19
16	Technical Data.....	19
17	Repair Return Form.....	20
11	Konformitätserklärung / Declaration of conformity.....	32

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1 User Instruktion

1.1 Important Instructions for your safety

- Every user must read and understand this manual completely before use. Only instructed users may operate the instrument. Failure to do so can result in serious injury or death.
- Follow general instructions for hazard prevention and general safety instructions, e.g. wear protection clothing, eye protection and gloves.
- This operating manual is part of the product. Thus, it must always be easily accessible.
- Enclose this operating manual when transferring the device to another place.
- If this manual is lost, please request another one. Please contact your dealer or

Ingenieurbüro CAT
M. Zipperer GmbH
Wettelbrunner Str. 6
D-79282 Ballrechten-Dottingen
Tel.: ++49-(0)7634-5056-800
Fax: ++49-(0)7634-5056-801
www.cat-ing.de
info@cat-ing.de

1.2 Danger symbols in this operating manual

The safety instructions in this manual appear with the following danger symbols and danger levels:

1.2.1 Danger symbols:

	Hazard point		Electrical shock
	Risk of fire		Explosion
	Crushing		Hot surface
	Bio hazard		Chemical hazard

1.2.2 Danger levels

	Will lead to severe injuries or death
	May lead to severe injuries or death
	May lead to light to moderate injuries
	May lead to material damage

1.2.3 Warning signs on the device

WARNING

This symbol indicates that it is imperative to read and understand the instruction manual prior to operating the instrument. Please highlight points which require special attention in your field of application so they are not overlooked. Disregard of warnings may result in impairment of serviceability as well as in physical harm to the user.

WARNING

This symbol indicates that there are **hot surfaces** on the device. It refers to the hotplate of the device which can reach up to 500°C during operation causing a hazard resulting in serious burns.

- Avoid contact with the hotplate, vessel and heating medium.
- Avoid storage of flammable substances near the instrument.

2 General safety warnings and instructions

	<p>⚠ DANGER Risk of explosion.</p> <ul style="list-style-type: none"> • Do not operate the device in the vicinity of highly flammable or explosive substances. The instrument is not explosion-proof. • Do not use this device for processing any substances which could generate an explosive atmosphere. • Do not use this device to process any explosive or highly reactive substances. • Do not use flammable or explosive substances near the instrument.
	<p>⚠ DANGER Electric shock as a result of penetration of liquid.</p> <ul style="list-style-type: none"> • Do not allow any liquids to penetrate the inside of the housing • Switch off the device and disconnect the power plug before starting cleaning or disinfection work. The On/Off Switch on the device does not disconnect the device from the power source. • Use only closed vessels. • Only plug the device back in if it is completely dry, both inside and outside.
	<p>⚠ WARNING Risk from incorrect supply voltage</p> <ul style="list-style-type: none"> • Only connect the device to an AC power source with a protective earth (PE). • Only connect the device to voltage sources which correspond to the electrical requirements on the type label.
	<p>⚠ WARNING Electric shock due to damage to device or mains cable</p> <ul style="list-style-type: none"> • Only connect the device to the mains supply if the device and the mains cable are undamaged • Only use devices that have been properly installed or repaired. • In case of danger, disconnect the device from the mains supply by pulling the power plug from the mains socket or by using the isolating device intended for this purpose (e.g. emergency stop switch)
	<p>⚠ WARNING Lethal voltage inside the device</p> <ul style="list-style-type: none"> • Do not open the device. • Ensure that the housing is always closed and undamaged so that no parts inside the housing can be contacted by accident. • The On/Off Switch on the device does not disconnect the device from the power source. Remove the plug from the AC power outlet to disconnect the instrument from the mains supply entirely. • Do not allow any liquids to penetrate the inside of the housing. • Repairs are only to be carried out by trained service technicians.
	<p>⚠ WARNING Damage to health due to corrosive or aggressive chemicals</p> <ul style="list-style-type: none"> • Observe all markings on the reagent bottles.

	<ul style="list-style-type: none"> ● Always check the instrument for leaks and air bubbles. Special attention should be directed to determine that all push-ons, threaded connections and suction tubes are firmly in place before beginning operation. Leaking solutions may endanger persons and materials ● Use proper connecting vessels, protective clothing and gloves. ● Avoid splashes ● When dispensing, maintain a physical distance between the instrument and the body. ● Dangerous and fuming chemicals must be dispensed in a fume hood. ● Only employ the instrument for the purpose intended by the manufacturer, and particularly within the resistance limits of the instrument. If in doubt, contact your supplier, or the manufacturer's factory representative at the phone number shown at the front page of this operating instruction. ● Always use the instrument in such a manner that neither the operator, nor any other person is endangered.
	<p>▲WARNING Damages to health due to infectious liquids and pathogenic germs.</p> <ul style="list-style-type: none"> ● When handling infectious liquids and pathogenic germs, observe the national regulations, the biological security level of your laboratory, the material safety data sheets and the manufacturer`s application notes. ● Wear personal protective equipment ● For comprehensive regulations about handling germs or biological material of the risk group II or higher, please refer to the “Laboratory Biosafety Manual” in its respectively current valid version from the World Health Organisation
	<p>▲WARNING Damage to health due to contaminated device and accessories In the following cases, sample material can be released:</p> <ul style="list-style-type: none"> - improperly closed vessels - unstable vessels - smashed glass vessels <ul style="list-style-type: none"> ● Only mix in closed vessels ● Observe the nationally prescribed safety environment when working with hazardous, toxic and pathogenic samples. Pay particular attention to personal protective equipment (gloves, clothing, goggles, etc.), extraction, and the safety class of the lab. ● Decontaminate the device and the accessories before storage and shipping.
	<p>▲WARNING Risk of fire</p> <ul style="list-style-type: none"> ● Do not use this device to process any highly flammable liquids

	<p>▲CAUTION Poor safety due to inadequate fixing of the unit</p> <ul style="list-style-type: none"> • Ensure that the unit is firmly attached to a solid stand.
	<p>▲CAUTION Poor safety due to incorrect accessories and spare parts.</p> <p>The use of accessories and spare parts other than recommended by Ingenieurbüro CAT, M. Zipperer GmbH may impair the safety, function and precision of the device. Ingenieurbüro CAT, M. Zipperer GmbH cannot be held liable or accept any liability for damage resulting from the use of incorrect or non-recommended accessories and spare parts, or from the improper use of such equipment.</p> <ul style="list-style-type: none"> • Only use accessories and spare parts recommended by Ingenieurbüro CAT, M. Zipperer GmbH
	<p>▲CAUTION Crush hazard due to moving parts</p> <ul style="list-style-type: none"> • Do not replace any consumables as long as the device is running.

3 Intended Use

The shaker is mainly to be used for mixing/shaking media in test tubes or micro plates. Always use suitable test tubes or micro plates.

The shakers are designed for use in chemical and biological laboratories of industrial enterprises, universities and pharmacies.

To ensure maximum service life, observe the specified ambient conditions (temperature and humidity) and ensure that the instrument is not exposed to a corrosive atmosphere.

▲WARNING It is the user's task to find out whether the device is suitable for the application. If in doubt clarify this with your dealer or directly with the manufacturer.

▲WARNING Please comply with all safety and accident-prevention regulations applicable to laboratory work.

4 General Information

The Shaker S 26 is designed in accordance with Safety Class 1 and built and tested in accordance with DIN EN 61010.

According to these regulations the unit is designed to meet the requirements for safe and correct operations. To maintain the proper safety and operational functions of the instrument the user should follow the instructions and safety guidelines in this manual.

5 Unpacking the Instrument

Unpack the instrument carefully and check to see that it is not damaged. It is important that any damage incurred in transport to be recognized at the time of unpacking. Notify your carrier or forwarding agent immediately in case of such damage.

5.1 Scope of delivery

Please check that the package contains the following:

Type	Description	Part No.
S 26	Shaker 230V Instruction Manual	60386-0000


Accessories to be ordered separately

Type	Description	Part No.
Cable	9-pole SUB-D-Cable	30275-0051
Adapter Cable	USB to RS232 (incl. Driver CD)	30244-0001
Shaking Attachment for Erlenmeyer flasks	6 x 250ml	60689-0010
	4 x 500ml	60689-0020
	2 x 1000ml	60689-0030

6 Set-up and Operation


6.1 The Front Panel Control

The keypad makes data entry easy and incorporates a 2-line alpha numeric LCD which shows all relevant data. The entry and changing of settings is made via the control knob in connection with the function keys and certain hot-keys.

Front panel view of the S 26

6.2 Description of the front panel functions

LCD-Display	<i>The LCD-Display shows all relevant system data.</i>
	Standby/On Button <i>Switch the device on or off.</i>
	Button: "Motor On/Off" <i>To switch on the motor this button needs to be pressed. Switching the motor off is also done by pressing this button.</i>
	Button: "Display RPM" <i>Shows the actual set motor speed.</i>
	Button: "Timer" <i>If required, pressing this key initiates the timer function which will shut down the S 26. The timer is calibrated to show: days: hours: minutes.</i>
	Speed Control Wheel: "Motor /RPM" <i>This wheel is used to set the speed for the shaking function. The set speed is displayed on the LCD display while this wheel is turned (opto-electronic speed control)</i>
	Control knob <i>This control knob is used to change different set values. Turning clockwise will increment values, while turning anti-clockwise will decrement.</i>
Status LED „Motor On/Off“	<i>This LED is illuminated when the motor is on. It indicates that the unit is trying to reach or keep the entered set speed.</i>
EDIT ARROW →:	<i>This appears in the display when data input is required .It stays on for three seconds, to allow the operator to set the values via the control knob. Should the arrow disappear during data entry, just press the key again.</i>

6.3 Set-up

Please put the unit on to an even surface. A minimum distance of 50 cm to inflammable materials should always be observed.

▲ DANGER

- The units are not to be used in environments with danger of explosion.

▲ WARNING

- Do not use this device to process any highly flammable liquids

▲ WARNING

- The unit is not to be used without supervision.

▲ WARNING

- When connecting the instrument to an AC power outlet, ensure that your local supply voltage matches that indicated on the instrument's rating plate.

6.3.1 Electrical Connection

▲WARNING Risk from incorrect supply voltage

- Only connect the device to an AC power source with a protective earth (PE).
- Only connect the device to voltage sources which correspond to the electrical requirements on the type label.

7 Operation Shaking

▲WARNING

Use only closed vessels before shaking to avoid any splashing out.

Before shaking, the user should take the following steps:

1. Put the required shaker attachment on the device and fix it with four star handle screws.
2. Put the required samples into the bracket. To prevent uneven vibrations, the vessels should be distributed uniformly along the shaker attachment.

Optimal shaking speed for water (H₂O):

Vessel	Filling amount	Shaking speed (1/min)
6 x 250 ml	each 180 ml	300 – 400
4 x 500 ml	each 300 ml	200 – 220
2 x 1000 ml	each 700 ml	190 – 250

To set the shaking speed and to switch the motor on, following steps are required:

▲WARNING

Make sure, that the device stands fixed on a base. Turn the speed control wheel slowly to the required shaking speed.

	<p>Turn the speed knob to set the desired shaking speed. Turn the wheel entirely to the left to switch the motor off.</p>
	<p>Press the Motor On/Off key to switch on the motor (Motor On/Off LED lights). To switch off the motor press the Motor On/Off key again (Motor On/Off LED goes out).</p>


7.1 Indexing mechanism

If the motor is switched off or speed control wheel is set to “0 rpm”, the motor automatically approaches the indexed position. As long as the unit is switched on but the motor is switched off, the motor stays at the indexed position.

8 Timer controlled operation

This device has an integrated clock timer. After an adjusted time, the timer switches the device into standby. If the timer expires, the device switches to standby mode (off).

The timer function can be switched on by the following way:

	<p>Press the Timer key.</p>
	<p>Turn the control knob to select the desired set value. The timer value is displayed in the format “dd:hh:mm”. For example the timer value “11d:22h:33m” means 11 days, 22 hours and 33 minutes. If the timer value is set to “00d:00h:00m” the timer function is disabled.</p>

Minimum shut off time: 1 Minute

Maximum shut off time: 99 days: 23 hours: 59 minutes

9 Error Messages

The S 26 has an internal self diagnostic system. This enables the device to detect errors to a large extent. Should one of the following errors occur the device will switch itself off. To enable the user to prevent the problem re-occurring, the S 26 will display any of the following error messages as it is switched on again.

Error message:	Probable cause of error:
TIMER expired	The device switched off automatically after the timer expired (this is not an error, just there to inform the user)

10 Serial Interface (RS232) of S 26

The serial interface (RS232) of the S 26 facilitates remote access of all functions as well as readout and change of all system parameters

A special daisy chaining feature facilitates chaining multiple units via their serial interfaces only through one PC-interface (RS232).

To setup a daisy-chain, the transmit line (TxD) of the PC is connected to the receive line (RxD) of the first controller. The transmit line of this controller is then connected to the receive line of the next controller in the chain. The transmit line of the last controller in the chain is returned to the receive line of the PC, which closes the link and forms the ring connection.

To address a specific controller in a daisy-chain, each controller carries a so called Slave Address, which can be any number from 1 to 255 (factory setting: 1). The address "0" is defined as the general call address on which all units in the chain will respond.


Parameters of the serial interface (1200,8,N,1):

Baudrate:	1200
Databits:	8 Bit
Parity:	none
Stopbits:	1

10.1 Connection between PC and S 26

The RS232 socket at the read side of the UNIT has the following pinout.

Pin	Description
2	transmit data TxD
3	receive data RxD
5	ground , GND

Socket of the serial interface (view on socket - 9 pin female Sub-D socket)

*For connection to a PC a standard 9pin prolongation cable is used (pins not crossed).

10.2 Format of an RS232 command

ADR , CMDCODE , PARAMETERLIST <CR>

Description:

ADR:	Slave Address of the desired unit
CMDCODE:	Command-code
PARAMETERLIST:	1 to 6 parameters separated by commas
CR:	The command string must be terminated by Carriage/Return (ASCII code 13)

10.3 RS232 - commands

Table 1:

CMD. CODE	Explanation	Parameter list	Range
WON	Set Status (ON/OFF control) of - Motor - Plate	1. Motor ON/OFF (0-> Off, 1-> On) 2. Plate ON/OFF (0-> Off, 1-> On)	0/1 0/1
RON	Read status of: - Motor - Plate	1. Dummy parameter to initiate transfer -> Controller sends in handshake: 1. Motor On/Off (0-> off, 1-> on) 2. Plate On/Off (0-> off, 1-> on)	1 0/1 0/1
RAC	Read actual	1. Dummy parameter to initiate transfer -> Controller sends in handshake 1. actual motor speed 2. actual plate temperature 3. actual probe temperature if connected, if not 'x' as return value)	1 0..1600 0..360 0..250, x
WSE	Write actual set values of: - Motor - Hotplate (plate) - external PT100 sensor (probe)	1. Setpoint of motorspeed in Rpm 2. Setpoint of hotplate temperature in °C 3. Setpoint of probe temperature in °C	0..1600 0..360 0..250
RSE	Read actual set values of: - Motor - Hotplate (plate) - external PT100 sensor (probe)	1. Dummy parameter to initiate transfer -> Controller sends in handshake: 1. Setpoint of motorspeed in Rpm 2. Setpoint of hotplate temperature 3. Setpoint of probe temperature	1 0..1600 0..360 0..250
WVO	Write/Set volume	1. Volume in ml	10..10000
RVO	Read volume	1. Dummy parameter to initiate transfer -> Controller sends in handshake: 1. programmed volume in ml	1 10..10000
WTR	Set / Write: -Timer - Ramp	1. Timer value in seconds. Set this value to 0 to disable the timer. 2. temperature ramp °C/min	0.. 2592000 1.. 450

CMD. CODE	Explanation	Parameter list	Range
	- Safety temperature	3. Safety temperature in °C	20.. 450
RTR	Read settings of: - Timer - Ramp - Safety temperature	1. Dummy parameter to initiate transfer -> Controller sends in handshake: 1. actual setting of the timer in seconds. If the timervalue is 0, the timer is disabled. 2. temperature ramp °C/min 3. Safety temperature in °C	1 0.. 2592000 1..450 20..450
WSM	Activate serial mode (→ enable/disable setpoint control via frontpanel)	1. Set this value to 0 to disable the setpoint control via RS232 commands, set this value to 1 to enable setpoint control via RS232 commands and disable the frontpanel control. (only used for units with motor)	0 / 1
WEE	Write parameters to non volatile EEPROM (e.g. slave-address, setpoints...)	1. Control parameter	2
WSA	Set slave address	1. New slave address	1.. 255

10.3.1 Format of the Controller Handshake

After receiving an RS232-command the unit will:

1. Send the received command to the next unit (or back to the PC if there is only one unit connected)
2. Answer with a handshake string, which is defined as follows:

ADR , "HS" , RETCODE , PARAMETERLIST CR

Before sending the next command, the host (PC) has to wait until the handshake string has been received from the unit(s).

Explanation:

ADR: Slave Address of the Controller sending the handshake

RETCODE: Error-code (see table 2 below)

PARAMETERLIST: 0 to 6 parameters (see table 1), each parameter is separated by a comma ", "

CR: the handshake as any command, is terminated by ASCII-code 13 (CR)

TABLE 2:

Return Code	Explanation	Parameter list
OK	command executed, no error	see table 1
UC	unknown command	None
PA	wrong parameter number (too few or too many parameters specified)	None
NA	command is not allowed in actual operation mode	actual operation mode
PR	at least one parameter is out of range	None
PL	at least one parameter is too long	None
DF	unknown data format	None

10.3.2 RS232-application example

Example 1: The S 26 should be programmed to operate with the following settings:

Shaking speed: 500 rpm

It is assumed, that the S 26 carries the RS232-slave address 1 = default

To program the unit with the above mentioned settings, the following commands have to be sent to the unit:

1,WSM,1 ; first enable set point control via RS232 (needs only to be send once, after power-up)
1,WSE,500,0,0 ; program set points (motor speed, hotplate, probe temp.).
1,WON,1,0 ; enable motor

The S 26 then sends the following strings to the next unit (or back to the PC):

1,WSM,1 ;The first string is the Echo of the received command.
1,HS,OK ;The second string indicates that the command was accepted and will be executed.

1,WSE,500,0,0 ;Echo of the second command.
1,HS,OK ;indicates that the command was accepted

1,WON,1,0 ;Echo of the third command.
1,HS,OK ;indicates that the command was accepted

Example 2: Reading the actual motor speed from the S 26

To readout the actual values from the unit the following command has to be sent:

1,RAC,1

The S 26 then sends the following two strings to the next unit (or to the PC):

1,RAC,1
1,HS,OK,550,0,x

The first string is the Echo of the received command.

The second string indicates that the command was accepted and will be executed it also contains the following information (see also table 1):

- the actual motor speed is 550 rpm

11 Maintenance and Cleaning

11.1 Cleaning

The outer casting is mostly 304 grade stainless steel together with anodised aluminium and a chemically resistant splash-proof membrane key pad. Therefore the unit can be easily cleaned with warm water and any appropriate liquid laboratory detergent.

11.2 Maintenance

In case of malfunction do not attempt to repair the unit. There are no user-serviceable parts in this instrument.

The shaker S 26 should only be opened and repaired by authorised service personnel. Any work on the electronics in the unit should only be carried out by knowledgeable, trained personnel. Any attempt by the user to repair the unit will immediately render the guarantee null and void. Please contact your local distributor in the event of a problem.

12 Dismantling, Transport and Storage

12.1 Dismantling

1. Switch the unit off.
2. Disconnect it from the mains supply.
3. Remove any glass beakers and any other equipment around the shaker.
4. Now you may remove the instrument from the working area.

12.2 Transport and Storage

Prior to transport:

Switch the instrument off and proceed with dismantling as described under “Dismantling”.

Place the instrument and its parts in its original packaging or another suitable container to protect it during transport. Close the packaging with adhesive tape.

Store the instrument in a dry environment. Please observe the specified conditions of the ambient (temperature and humidity).

Do not subject the instrument to mechanical shocks or vibration during transporting it.

In case you do not use the original packaging please mark the box with the following notes:

- Glass symbol (handle with care, fragile)
- Umbrella (keep dry)
- Content (list of content)
- Storage ambient:
Max. ambient temperature : RT to +40°C
Max. humidity: 80%

13 Disposal

Please dispose of used instruments and defective components at your local recycling collection point. Prior to disposal, sort according to materials: Metal, glass, plastic, etc. Also be sure to dispose of the packing material in an environmental-friendly manner.

14 Warranty and Liability

The manufacturer agrees to either repair, or replace, at the manufacturer's discretion, any defects in materials or workmanship which develop within 24 months of the delivery of this product to the original user. In the event of replacement, the replacement unit will be guaranteed for the remainder of the original twenty-four (24) months period or ninety (90) days, whichever is longer.

If this product should require service, contact your local distributor or manufacturer for necessary instructions.

This guarantee will not apply if the defect or malfunction was caused by accident, neglect, unreasonable use or fitness for a particular purpose, which extend beyond the description and period set forth herein.

The manufacturer's sole obligation under this guarantee is limited to the repair or replacement of a defective product and the manufacturer shall not, in any event, be liable for any incidental or consequential damages of any kind, resulting from use or possession of the product.

▲WARNING

The user has to determine, if the instrument is suitable for his specific application. If there are any further queries, contact your local dealer or the manufacturer direct.

15 Technical Data

Type	S 26
Movement type	horizontal rotatory
	Orbit 3 mm
Display	backlit LCD display
Motor	capacitor motor
Motor power (output)	4 W
Motor control	PID (Microprocessor)
Motor speed	60 - 1100 1/min
Serial Interface	RS232
additional functions	Timer, optional Indexing Mechanism
Supply voltage	230V/50 Hz
Power	18 W
Duty cycle	100%
Allowable ambient temperature for operation	5 - 40°C
Humidity	80 % relative humidity
Safety class (DIN EN 60529)	IP42
Dimensions (W x H x D) [mm]	275 x 245 x 150 mm
Weight [kg]	6,9 kg
Order No.	60386-0000
Type	S 26

16 Repairs

▲WARNING

- When returning instruments for repair that have come into contact with hazardous substances, please provide precise information on the relevant medium
- Take protective measures to ensure the safety of our receiving and maintenance personnel
- Mark the package as appropriate for hazardous materials.

Add the completed Repair return form to the unit when sending back for repair.

17 Repair Return Form

CONTACT/USER INFORMATION

Contact:	Phone No.
Fax No.	Email:

Billing:	Shipping:
Company	Company
Address	Address

INSTRUMENT INFORMATION

Model	Serial No.	
Please describe all problems/malfunctions		
Operating Conditions (please fill in if applicable)		
Ambient Temp.	Humidity	Speed
Load	Volume	Viscosity
Temperature in °C	Sample Temperature	Operating Time
Sample Description*		

*NOTE: If the instrument was exposed to hazardous material, it must be decontaminated BEFORE returning it to Ingenieurbüro CAT, M. Zipperer GmbH and an MSDS for hazardous material must be included with the instrument.

RETURN SHIPPING

UPS <input type="checkbox"/>	Air Parcel Post <input type="checkbox"/>	Collect** <input type="checkbox"/>	Other** <input type="checkbox"/>
------------------------------	--	------------------------------------	----------------------------------

**Your account number is required for UPS collect respectively the address and contact of your preferred forwarder if you choose any other transport means.

PACKAGING INSTRUCTIONS TO RETURN AN INSTRUMENT FOR REPAIR
--

- ✓ Remove all accessories (e.g. homogenizer tools, stirring paddles) from the instrument
- ✓ Clean excess testing material off the instrument/accessory
- ✓ Include MSDS sheets for all hazardous materials used with this instrument
- ✓ Pack the instrument in its original box. If the box is not available, take care to wrap the instrument and accessories with enough material to support them.
- ✓ DO NOT send pedestal stand unless there is a problem with the upright rod, clamp or base. If there is a problem with the stand remove the upright rod from the base and individually wrap each item to avoid contact with the instrument. (Applicable for overhead stirrers and homogenizers.)
- ✓ Pack the instrument and related items in a strong box for shipping. Mark the outside of the box with handling instructions.

Example: "Handle with care" or "Fragile- Delicate Instrument" and send to:

Ingenieurbüro CAT, M. Zipperer GmbH, Wettelbrunner Str. 6, D-79282 Ballrechten-Dottingen

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We declare under our sole responsibility that this product (see table) corresponds to the regulations and conforms with the standards or standardized documents:

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The agreement with further valid guidelines/regulations following for the product is explained:

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Applied one harmonized standards in particular:

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In the case of a modification of the unit which has not been agreed on with us, this declaration becomes null and void and the warranty expires.

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